

**24-Plex Carboxyl Barcoded Magnetic Beads Protein Coupling Kit**  
**Protein/Antibody Coupling and Immunoassay Instructions**

**Part Number: 64-R0112-PIN**

For Exclusive Use with  
Applied BioCode 1000A or 2000 Analyzer

---

***For Research Use Only - Not for Use in Diagnostic Procedures***

## **INTRODUCTIONS**

### **Intended Use**

The 24-plex Carboxyl Barcoded Magnetic Bead Protein Coupling kit is a research use only kit for the purpose of enabling users to create their own assays using their own antibodies or proteins. The kit contains components required to build assays for two 96-well plates. It also lists other reagents and tools that are required to complete assay development.

### **Principle**

The Applied BioCode System is a flexible multiplexing platform for detecting and analyzing targets using the Barcoded Magnetic Beads (BMB). A wide variety of assay types, such as nucleic acid hybridization assays, and immunoassays are performed in an aqueous, both quickly and efficiently.

The BioCode Analyzer and BMB technology offers multiplex capability for simultaneous detection up to 128 different analytes within a single sample.

With BMB technology, nucleotide/antibody reactions take place on the surface of BMB. For each sample, target specific capture nucleotide/antibody probes are covalently linked on to a specific set of BMB: (1) biotin labeled targets are captured by the BMB-bound oligo probes in a hybridization assay, (2) analyte/antigen captured by a specific capture probe coupled on BMB will in turn immobilize a biotin labeled detection antibody in a solid-phase immunoassay. Finally, the Streptavidin R-phycoerythrin conjugate is added to the samples for quantitative fluorescence detection.

The BioCode software completes the data collection and reports the results in a matter of minutes.

## STORAGE and STABILITY

1. Kit is shipped at ambient temperature. Store remaining kit components at room temperature (RT, 15 to 30°C).
2. All components are guaranteed up to the expiration date found on the packaging label if handled properly.

## SAFETY

1. Please refer to the product's Material Data Safety Sheet (MSDS) for safety information concerning this product.
2. Avoid reagent contact with skin. If contact is made, thoroughly wash the area with water.
3. Handle all specimens in accordance with Universal Precautions. (CDC, Update: Universal Precautions for Prevention of Transmission of Human Immunodeficiency Virus, Hepatitis B Virus, and other Bloodborne Pathogens in Health-Care Settings. MMWR 1988; 37:377-388. OSHA standards: 1910.1030 (d) (l) and 1910.1030 (d) (3).)
4. Waste must be classified and disposed of in accordance with all Federal, State, and Local environmental regulations.

## NECESSARY REAGENTS, MATERIALS and EQUIPMENT

### Kit Contents

Reagent Description	Quantity
P-Carboxyl Barcoded Magnetic Beads (BMB), 24-plex. The BMB's surface is functionalized for <u>protein/antibody</u> coupling.	24 tubes, ~10,000 beads/tube in 0.5 mL PBST

Part #	BMB #	Included (check mark)	Part #	BMB #	Included (check mark)
64-R0112-a	0 to 23		64-R0112-d	72 to 95	
64-R0112-b	24 to 47		64-R0112-e	96 to 119	
64-R0112-c	48 to 71				

<b>Activation Buffer</b> , 15% Ethanol in 100 mM MES, 0.05% ProClin, pH 6.0	1 vial, 20 mL
<b>Coupling Buffer</b> , 140 mM Guanidine-HCl, 100 mM MES, 0.05% ProClin 950, pH 5.0	1 vial, 20 mL
<b>PBS</b> , 1X, 0.05% ProClin 950, pH 7.4	1 vial, 20 mL

<b>PBST</b> , 0.05% Tween-20, in PBS w/ 0.05% ProClin 950, pH 7.4	1 bottle, 450 mL
<b>TRIS Buffer</b> , 50 mM TRIS-HCl, 0.05% ProClin 950, pH 7.4	1 vial, 20 mL
<b>Tween/NaCl</b> , 1 M NaCl, 0.01% Tween-20, 0.05% ProClin, pH 7.5	1 vial, 20 mL
<b>Detection Buffer</b>	1 bottle, 60 mL
96-well Pate/Cover	2 plates and 1 cover

**Reagent Preparation** (prepares fresh just before use, discard after use)

<b>Reagent Description</b>	<b>Protocol</b> (sufficient volume for 1 pate)
<b>Blocking Buffer</b> , 1% BSA, 1% NFDM in PBS	Dissolve 60 mg BSA and 60 mg NFDM in 6 mL PBS
<b>Assay Buffer</b> , 1% BSA in PBST	Dissolve 40 mg BSA in 40 mL PBST
<b>EDC</b> , 50 mg/mL	Add 200 $\mu$ L Activation Buffer to EDC vial
<b>Sulfo-NHS</b> , 50 mg/mL	Dissolve 25 mg Sulfo-NHS in 0.5 mL Activation Buffer
<b>SA-PE</b> , 2 $\mu$ g/mL	Dilute 10 $\mu$ L of 1 mg/mL SA-PE in 5 mL Assay Buffer

**Reagent Components and Supplies Not Provided**

<b>Description</b>	<b>Reference</b>
Sulfo-NHS (N-hydroxysulfosuccinimide)	ProteoChem, catalog # C1102 or equivalent
EDC, 10 mg	ProteoChem, catalog # C1100 or equivalent
BSA (Bovine Serum Albumin), protease, DNase free powder	Equitech Bio, catalog # BAH67 or equivalent
NFDM (non-fat dry milk)	Village Farm, Strum Foods, Inc or equivalent
SA-PE (Streptavidin R-phycoerythrin) with BSA, 1 mg/mL	MOSS Inc., catalog # SAPE-001 or equivalent
EZ-Link Sulfo-NHS Biotinylation Kit	Thermo Scientific, catalog # 21435
Zeba Spin Desalting Columns	Thermo Scientific, catalog # 89882

**Required Equipment**

<b>Description</b>	<b>Reference</b>
BMB Analyzer	Applied BioCode 1000A or 2000
Magnetic Stand, 12-positon	Promega, catalog # Z5342
Magnetic 96-well Plate Separator	Applied BioCode, Part # 01-M0001
Micro-centrifuge, bench-top	Galaxy Mini Star, VWR or equivalent
Vortex Mixer	VWR
Thermo Shaker, 96-well plate	Vortemp 56, Labnet or equivalent


Shaker, tube	BioShaker XP, Q Instruments or equivalent
Pipettor: P-20, P-200, P-1000, Multi-channel	

## RECOMMENDED PROTOCOLS

### Antibody Coupling

#### Technical Notes

- The process is for coupling antibodies or similar proteins to carboxylated BMB.
  - The optimal binding capacity on the BMB may depend on the robustness of the antibody and the pH of Activation and Coupling Buffers. This kit includes the Activation Buffer (pH 6) and Coupling Buffer (pH 5) which is suitable for most antibodies.
  - If the antibody is stored in buffers containing free amines, such as TRIS or Glycine, should be removed using a desalting column before coupling.
  - For coupling, we recommend using 2 to 10  $\mu\text{g}$  of antibody/protein per 10,000 beads.
  - Prepare the EDC solution immediately before use.
1. Spin down and place the BMB tubes into the Magnetic Stand for 1 to 2 minutes.  
Note: Coupling up to 12 reactions can be performed simultaneously.
  2. Remove the supernatant using a pipette.
  3. Wash the BMB twice - add 200  $\mu\text{L}$  of Activation Buffer into the BMB tube, gently vortex the tube for 3 to 5 seconds. Repeat Steps 1 and 2.
  4. Repeat Step 3 once for a total of two washes.
  5. Add 76  $\mu\text{L}$  of Activation Buffer into the BMB tube.



6. Add 14  $\mu\text{L}$  of 50 mg/mL Sulfo-NHS into the BMB tube, gently vortex the tube for 3 to 5 seconds.
7. Add 10  $\mu\text{L}$  of 50 mg/mL EDC into the BMB tube, gently vortex the tube for 3 to 5 seconds; centrifuge the BMB tube for 8 to 10 seconds.
8. Place the BMB tube in a shaker for 20 to 30 minutes, mix at 1200 to 1500 rpm, at RT.



9. Remove the Activation Buffer from the BMB tube, follow Steps 1 and 2.
10. Wash the BMB twice – add 200  $\mu\text{L}$  of Coupling Buffer into the BMB tube, gently vortex the tube for 3 to 5 seconds. Repeat Steps 1 and 2.
11. Repeat Step 10 once for a total of two washes.
12. Add 100  $\mu\text{L}$  of Antibody (or Protein) diluted in Coupling Buffer into the BMB tube.
13. Place the BMB tube in a shaker for 2 hours, mix at 1200 to 1500 rpm, at RT.
14. Remove the supernatant, follow Steps 1 and 2.
15. Wash the BMB once – add 200  $\mu\text{L}$  of Tween/NaCl into the BMB tube, gently vortex the tube for 3 to 5 seconds. Repeat Steps 1 and 2.
16. Add 200  $\mu\text{L}$  of PBST into the BMB tube.
17. Place the BMB tube in a shaker for 13 to 15 minutes, mix at 1200 to 1500 rpm, at RT.
18. Remove the supernatant, follow Steps 1 and 2.
19. Wash the antibody/BMB tube twice – add 200  $\mu\text{L}$  of TRIS Buffer into the BMB tube.
20. Place the BMB tube in a shaker for 30 minutes, mix at 1200 to 1500 rpm, at RT.
21. Remove the supernatant, follow Steps 1 and 2.
22. Repeat Steps 19 to 21 except shake the BMB tube for 10 minutes. For a total of two washes.


23. Wash the antibody/BMB tube once – add 200  $\mu$ L Tween/NaCl into the BMB tube.
24. Place the BMB tube in a shaker for 3 to 5 minutes, mix at 1200 to 1500 rpm, at RT
25. Remove the supernatant, follow Steps 1 and 2.
26. Wash the antibody/BMB tube once – add 200  $\mu$ L of Blocking Buffer, gently vortex the tube for 3 to 5 seconds. Repeat Steps 1 and 2.
27. Wash the antibody/BMB tube once – add 200  $\mu$ L of Blocking Buffer into the BMB tube.
28. Place the BMB tube in a shaker for 30-60 minutes, mix at 1200 to 1500 rpm, at RT.
29. Remove the supernatant, follow Steps 1 and 2.
30. Wash the antibody/BMB tube twice - add 200  $\mu$ L of Assay Buffer into the BMB tube, gently vortex the tube for 3 to 5 seconds. Repeat Steps 1 and 2.
31. Repeat Step 30 once for a total of two washes.
32. Add 650  $\mu$ L of Assay Buffer into the BMB tube.
33. Store BMB tube at 2 to 8°C until needed.

### **Biotinylation of Detection Antibodies**

1. Recommended to use the Thermo Scientific Sulfo-LC Biotinylation kit.
2. Store the Biotin-Antibody in Assay Buffer at 2 to 8°C.
3. The Biotin-Antibody working concentration should be optimized.

### **Immunoassay Protocol**

1. Recommend to use approximately 50 beads/plex per well.
2. Calculate the BMB volume needed for number of run(s) with 50 beads/plex per well.
3. Combine the BMB volume of each plex (analyte) into a master-mix tube.
4. As the example above, vortex each Stock Antibody/BMB tube until the beads are suspended, immediately transfer 0.31 mL into a 10-mL tube (Master-mix).  
Continues this transfer process until the desired plex numbers are in the Master-mix tube.
5. If necessary, add Assay Buffer to the Master-mix tube to a final volume approximately 7.7 mL.




Notes:

- The 7.7 mL final volume is for running one full plate.
- Calculate the make-up volume if it is less than a full plate. Example if only running 30 wells assay:

$$\begin{aligned} \text{Final Vol.} &= 75 \mu\text{L} \times \text{number of wells} \times 1.05 \\ &= 75 \mu\text{L} \times 30 \times 1.05 = \mathbf{2.4 \text{ mL}} \end{aligned}$$

6. Transfer 75  $\mu\text{L}$  of Antibody/BMB from the Master-mix tube into corresponding wells.
7. Place the BMB plate on the Magnetic Plate Separator for 1 to 2 minutes.
8. Remove the supernatant using a multi-channel pipette.



9. Add 50  $\mu\text{L}$  of Calibrators / Unknown Samples diluted in Assay Buffer into the corresponding wells.
10. Place the covered BMB plate in a shaker for 1.5 – 2 hours, mix at 900 to 1000 rpm, at RT.



11. Remove the supernatant from the BMB plate, follow Steps 7 and 8.
12. Wash the BMB twice – add 200  $\mu\text{L}$  of PBST per well, follow Steps 5 and 6.
13. Add 50  $\mu\text{L}$  of Biotin-Antibody diluted in Assay Buffer into each well.
14. Place the BMB plate in a shaker for 1 hour, mix at 900 to 1000 rpm, at RT.
15. Remove the supernatant from the BMB plate, follow Steps 5 and 6.
16. Wash the BMB twice as Step 10.
17. Add 50  $\mu\text{L}$  of 5  $\mu\text{g}/\text{mL}$  SA-PE into each well,
18. Place the covered BMB plate in a shaker for 30 minutes, mix at 900 to 1000 rpm, at RT.
19. Remove the supernatant from the BMB plate, follow Steps 7 and 8.
20. Wash the BMB 3 times as Step 10.
21. Add 200  $\mu\text{L}$  of Detection Buffer into each well.
22. Place the covered BMB plate in a shaker for 2 minutes, mix at 700 rpm, at RT.
23. Measure the fluorescence intensity in the BioCode 1000A or 2000 Analyzer.

## TECHNICAL SERVICE AND ORDERING INFORMATION

1. Customer Service: 562-801-0050
2. Ordering: Please call 562-801-0050 ext. 253 or email your orders to [orders@ApBioCode.com](mailto:orders@ApBioCode.com)
3. Bulk reagent order:



<b>Description</b>	<b>Size</b>	<b>Part Number</b>
Activation Buffer, pH 6.0	100 or 250 mL	44-A0501
Coupling Buffer, pH 5.0	100 or 250 mL	44-C0501
PBS, 1X, pH 7.4	450 or 950 mL	44-P0501
PBST, pH 7.4	450 or 950 mL	44-P0502
TRIS Buffer, pH 7.4	100 or 250 mL	44-T0502
Tween/NaCl, pH 7.5	100 or 250 mL	44-T0503
Detection Buffer	500 or 1000 mL	44-D0002
96-well Plate	10 plates/pack, 2 plate covers	01-P0002, 01-P0004
P-Carboxyl BMB	50,000 beads/tube	44-B0112

**DISTRIBUTED BY:**

**Applied BioCode, Inc.**  
10020 Pioneer Blvd., Suite 102  
Santa Fe Springs, CA 90670 USA  
Tel: (562) 801-0050  
[www.ApBioCode.com](http://www.ApBioCode.com)