

MULTIPLEX GASTROINTESTINAL PATHOGEN DETECTION WITH A USER FRIENDLY, HIGH-THROUGHPUT SYSTEM

Bansari Shah, Melissa Henrie, Jakob Kirchner, Deepali Shinde, Tad Kawashima, Michael Aye*
Applied BioCode, Inc., Santa Fe Springs, CA

Revised Abstract

Gastroenteritis is the second most common illness after the common cold. Globally, diarrhea accounts for approximately 2 million deaths in children under 5 years old each year, or 19% of total child deaths.^{1,2} High-throughput multiplex assays are desirable for rapid identification of the pathogens that can cause outbreaks of diarrhea and for infection control in healthcare settings.

Using the proprietary barcoded magnetic bead (BMB) technology, Applied BioCode has developed a prototype molecular diagnostic assay for detection of gastrointestinal (GI) pathogens including bacteria, viruses, and parasites. The GI pathogen assay is highly specific and did not show cross-reactivity with organisms tested in the study. Analytically, the assay can detect 10 copies of DNA targets per PCR reaction, and preliminary limit of detection was as low as 10³ CFU/mL for *C. difficile* and 10² PFU/mL for adenovirus 41. Evaluation of 43 clinical specimens previously tested by validated real-time PCR assays showed 100% concordance for *C. difficile* toxin and norovirus. In parallel with the multiplex assay development, we are also developing a prototype automated system in a 96-well format. This system will integrate PCR with post-PCR processing and detection, thereby simplifying the workflow and reducing hands-on time for the users.

In this study, we tested 234 stool specimens with the BioCode GI assay and demonstrated comparable results to a FDA-cleared multiplex test performed in a clinical laboratory.

Methods

Clinical specimens were subjected to vortexing with beads in the easyMAG lysis buffer prior to automated extraction with NucliSENS easyMAG® (bioMerieux).

Following extraction of nucleic acids from stool specimens, DNA and RNA targets were amplified in a one-step RT-PCR in either 9700 or 7500 Fast thermal cycler (Applied Biosystems). PCR products were captured by target-specific probes that were coupled to BMBs with unique barcodes. The presence of captured target sequence on an individual BMB was detected by streptavidin-phycoerythrin (SA-PE) conjugate. Qualitative results for each target were determined by a median fluorescence index (MFI) signal relative to assay cutoffs.

Target specific oligonucleotide probes were synthesized with a 5'-amino modification and a linker. Each probe was conjugated to BMBs with a unique barcode via amide linkage. Probes coupled to BMBs were combined into assay specific probe mix used for capturing biotin-labeled target DNA. Target capture (hybridization) and SA-PE labeling reactions were performed in a Thermo Shaker (PlexBio Co.). Unbound molecules were removed from BMB by washing with a ELX50 magnetic plate washer (BioTek). BMB imaging and fluorescence signal detection was performed on BioCode 1000A reader (Applied BioCode, Inc.).




Figure 1. Schematic presentation of a biotinylated target DNA captured by a target-specific probe (red line) coupled to BMB (gray rectangle), and labeled by SA-PE (*S)

BioCode GI Panel Targets

Table 1. Organisms and toxins targeted by BioCode GI Panel

Bacteria	
• <i>Clostridium difficile</i> (toxin A and B)	
• <i>Campylobacter</i> spp.	
• <i>Salmonella</i> spp.	
• <i>Shigella</i> spp./ Enteroinvasive <i>E. coli</i>	
• <i>E. coli</i> O157	
• Enterotoxigenic <i>E. coli</i> (LT/ST)	
• Shiga toxin producing <i>E. coli</i> (stx1/stx2)	
• Enterococcal <i>E. coli</i>	
• <i>Aeromonas</i> spp.	
• <i>Vibrio</i> spp. (<i>V. parahaemolyticus</i>)	
Viruses	
• Norovirus GI/GII	
• Adenovirus 40/41	
• Rotavirus A	
Parasites	
• <i>Cryptosporidium</i> spp.	
• <i>Giardia intestinalis</i>	
• <i>Entamoeba histolytica</i>	
Internal Control (MS2)	

Barcoded Magnetic Beads (BMB)




Figure 2. BMBs extracted from wafers are activated in suspension, coupled to probes, and utilized for target capture in microtiter plates.

BioCode-3000 HTS





Figure 3. Workflow for BioCode assays includes automated extraction and PCR set up prior to the steps automated by BioCode-3000 HTS (Top). Computer-aided drafting (CAD) image of the prototype instrument (Bottom).



Results




Figure 4. Sample data from the BioCode assay with analytical targets. MFI values (bars) for each target (X-axis) are plotted for each probe.

Table 2. Microorganisms tested for cross reactivity with BioCode assay

Organisms tested for Cross Reactivity		
<i>Acinetobacter baumannii</i>	<i>Acinetobacter lwoffii</i>	<i>Alcaligenes faecalis</i>
<i>Bacillus cereus</i>	<i>Citrobacter freundii</i>	<i>Enterobacter aerogenes</i>
<i>Enterobacter cloacae</i>	<i>Enterococcus faecalis</i> (vanB)	<i>Enterococcus faecium</i> (vanA)
<i>Escherichia coli</i>	<i>Klebsiella oxytoca</i>	<i>Klebsiella pneumoniae</i>
<i>Lactobacillus acidophilus</i>	<i>Listeria monocytogenes</i>	<i>Proteus mirabilis</i>
<i>Providencia stuartii</i>	<i>Pseudomonas aeruginosa</i>	<i>Serratia marcescens</i>
<i>Staphylococcus aureus</i>	<i>Staphylococcus epidermidis</i>	<i>Stenotrophomonas maltophilia</i>
<i>Streptococcus galactiae</i>	<i>Vibrio cholerae</i>	<i>Yersinia enterocolitica</i>

No cross-reactivity was observed for the BioCode GI assay with the microorganisms listed above.

Table 3. Preliminary Limit of Detection (LoD) for the BioCode GI assay

Target	Preliminary Limit of Detection	
	Preliminary LoD	Average MFI
Adenovirus 41	1x10 ² PFU/mL	4898
Norovirus GI (recombinant)	7x10 ³ TCID ₅₀ /mL	6937
Norovirus GII (recombinant)	7x10 ³ TCID ₅₀ /mL	5672
<i>C. difficile</i>	1x10 ³ CFU/mL	2624
<i>Salmonella enterica</i> (typhimurium)	7x10 ³ CFU/mL	4256
<i>Shigella flexneri</i>	7x10 ³ CFU/mL	5721
<i>E. coli</i> O157	7x10 ³ CFU/mL	6130
<i>Aeromonas hydrophila</i>	7x10 ³ CFU/mL	7103
<i>Vibrio parahaemolyticus</i>	7x10 ³ CFU/mL	7443
<i>Cryptosporidium</i>	7x10 ³ cells/mL	6736
<i>Giardia intestinalis</i>	1x10 ² cells/mL	1241

Table 4. Comparison of the BioCode assay with validated real-time PCR for *C. difficile* and noroviruses (GI & GII).

BioCode GI Panel	Real-Time PCR		
	Pos	Neg	Total
	Total	10	33

BioCode GI Panel	Real-Time PCR		
	Pos	Neg	Total
	Total	8	35

Results of the BioCode assay showed excellent correlation with validated real-time PCR assays for *C. difficile* and noroviruses.

Results

Table 5. Comparison of the BioCode GI assay with Luminex GPP assay (n=234 stool specimens)

Target Pathogen	Positive results reported by	
	BioCode GI	Luminex GPP
<i>Clostridium difficile</i>	67	65
<i>Salmonella</i> spp.	5	7
<i>Shigella</i> spp./ EIEC	3	4
ETEC	0	0
STEC	0	0
<i>E. coli</i> O157	1	1
EAEC	5	N/A
<i>Campylobacter</i> spp.	1	1
<i>Aeromonas</i> spp.	1	N/A
<i>Vibrio parahaemolyticus</i>	0	N/A
Norovirus (GI & GII)	11	10
Rotavirus A	10	9
Adenovirus F	5	N/A
<i>Cryptosporidium</i> spp.	1	1
<i>Giardia intestinalis</i>	0	0
<i>Entamoeba histolytica</i>	0	N/A
Total positive results	110	98
Invalid results (initial run)	1 (8)	1 (27)*

*Repeat testing for Luminex assay was performed with 1:10 dilution of the sample.

Co-infections: The BioCode assay reported 10 out of 234 samples (4.3%) as positive for more than one pathogen.

Table 6. Comparison of BioCode assay with Luminex GPP for *C. difficile*.

BioCode GI Panel	Luminex GPP Results		
	Positive	Negative	Total
	Total	59	161

Positive Agreement		90.8%
Negative Agreement		95.3%

Summary & Conclusions

- BioCode GI assay is shown to detect intended targets including bacteria, viruses and parasites.
- BioCode GI assay does not detect other organisms tested for cross-reactivity.
- Preliminary Limit of Detection (LoD) of the assay was determined at 10³ CFU/mL for *C. difficile* 10² PFU/mL for adenovirus 41 and 10² cells/mL for Giardia.
- Clinical performance with 234 stool specimens showed comparable results to the Luminex xTAG® GPP assay.
- BioCode-3000 HTS is designed for batch testing of multiplex assays in 96-well format with limited hands-on time and reduced contamination risk.

Acknowledgements

The authors thank Diatherix Laboratories, Inc. for providing some of the itered organism stocks used for preliminary LoD of the BioCode GI assay.

References

- Boschi-Pinto C, Velebit L, Shubuya K. Estimating child mortality due to diarrhea in developing countries. Bull World Health Organ. 2008; 86(9):710-7.
- Clark B, McKendrick M. A review of viral gastroenteritis. Curr Opin Infect Dis. 2004; 17(5):461-9.